

HPLC Column Heater

FS-10,20

FS-01

Widely used in gravimeter, microanalysis, colloid dispersion and bacteria-free experiments, this HPLC filter holder can filtrate and degas the mobile phase of HPLC analysis. Besides playing a critical role in lifetime extension of apparatus and HPLC column, the filter holder can enhance measurement accuracy.

The filtration bottles have high pressure endurance, fabulous interchange ability, and well-proportional wall. They are made according to international standard sizes so that they can match most products of foreign brands. In addition, the newly designed solvent pickup adaptor can filtrate or collect solvent directly from solvent bottle through film of filter. This new design makes the process of filtrating become easier and it can dramatically increase effectiveness and reduce air pollution.

Features

- High-quality and hard-glass materials that can endure more than 200 temperature difference.
- Elegantly made with well-proportional wall.
- Good endurance of high pressure and fabulous interchange ability.
- Made according to international standard sizes that can match most foreign brand products.

Packing List

FS-10, 20	FS-01
Glass Funnel: 47mm, 300ml	Glass Funnel: 47mm, 300ml
Plug	Silicone Plug
Fritted Glass Support Base: 47mm	Aluminum Clamp: 47mm
Aluminum Clamp: 47mm	Fritted Glass Support Base: 47mm
Erlenmeyer Flask: 1000ml (FS-10), 2000ml (FS-20)	Glass Cover
	Beaker
	Teflon Base

CROMTECH INDIA

1st Floor Begum Niwas, Opp Mercantile, Bank Society Colony, 5th Natwar Nagar Road, Jogeshwari (East), Mumbai 400 060.

Maharashtra - India.

Mobile: +91 98693 82134 / 81088 13320

Telefax: +91 22 2838 6028

Email: cromamumbai@mail.com / cromtechindia@mail.com / info@cromtechindia.com / sales@cromtechindia.com

Website: www.cromtechindia.com